


## • Ease of use

Handle adopting four-link structure with clear stops, achieve easier sampling and comfortable operation





## • Enlightened jaw

Stainless steel stamping head, laser welding, more reliable connection







## Single-use



No need to wait for reprocessing, Always available for your needs

## Achieve excellent passability



The tube wrapped in ultra-smoot plastic, can reduce the damage to the working channel

Order Number	Jaw Outer Diameter (mm)	Working Channel (mm)	Working Length (mm)	Jaw Type	Needle	Coating	Application	Jaw Type
FB-12SX-B1	1.0	≥1.2	700	Tooth Alligator Jaw	No	No	G,D	
FB-12C-B1	1.0	≥1.2	1050	Tooth Alligator Jaw	No	No	A 0	
FB-12D-B1	1.0	≥1.2	1100	Tooth Alligator Jaw	No	No	O <sup>H</sup> O	
FB-12E-B1	1.0	≥1.2	1200	Tooth Alligator Jaw	No	No	010	
FB-17C-B1	1.5	≥1.7	1050	Tooth Alligator Jaw	No	No	M ©	
FB-20C-A1	1.8	≥2.0	1050	Fenestrated elongated cup	No	No	A ©	=#
FB-20C-A3	1.8	≥2.0	1050	Fenestrated elongated cup	No	Plastic coating	A 0	
FB-20C-B1	1.8	≥2.0	1050	Tooth Alligator Jaw	No	No	Ø\ <b>©</b>	
FB-20C-B3	1.8	≥2.0	1050	Tooth Alligator Jaw	No	Plastic coating	M ©	

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

