

THE RIGHT MEASUREMENT, RIGHT AWAY

See your confidence in a single, automated IOP measurement with ActiView™ Positioning System

Tono-Vera[™] Vet Tonometer







Measure eye pressure effortlessly with Tono-Vera[™] Vet

The all-new Reichert® Tono-Vera™ Vet Tonometer featuring the advanced, patent-pending ActiView™ Positioning System, quickly guides you to the apex of the cornea, giving you confidence in your intraocular pressure (IOP) readings. Tono-Vera Vet automatically measures when aligned, providing a more objective and repeatable result in just one measurement.

Designed with the veterinary professional in mind, the innovative Tono-Vera Vet system, which includes a base that stores and dispenses $Ocu-Dot^{\mathbb{M}}$ Tonometer Probes, is a perfect solution for every practice. Tono-Vera Vet is the only veterinary tonometer with Bluetooth*, making transferring data easier than ever.

Reichert® is Technology, Tradition, and Trust







A more objective measurement with ActiView[™]

ActiView $^{\mathbb{M}}$ Positioning System features a full color view of the eye combined with intuitive, interactive alignment prompts to guide the user precisely to the center of the cornea every time. When proper alignment is achieved, Tono-Vera $^{\mathbb{M}}$ Vet automatically measures, providing you with a result in just a single measurement. Intuitive color-coded rings around the final IOP result indicate the reliability of the measurement.

All the important information, always on display

Tono-Vera Vet has a comprehensive yet simple user interface

- 1 Selected Species
- 2 Battery Status
- 3 Bluetooth® Connectivity
- 4 ActiView Target
- 5 Probe Angle Indicator
- 6 Number of Measurements
- 7 Measure Mode
- 8 OD (Right Eye)
- 9 Access to Menu
- 10 OS (Left Eye)
- 11 Select OD/OS Buttons





Tono-Vera[™] Vet is as innovative as it is thoughtfully designed

Lightweight, Ergonomic, and Calibration-free

The perfectly placed operator-facing display and lightweight, symmetrical design allow for comfortable left or right-handed operation. Intuitive button placement provides easy access to all settings and measurement functions. Always ready to measure with no routine calibration required.

Rechargeable, Standard AA, or Both

Choosing Tono-VeraTM Vet Rechargeable Model means fewer batteries in landfills and a device that is always ready to use. Tono-Vera Vet AA Model goes anywhere and allows for a quick and easy battery change with no additional tools required. Rechargeable and AA battery packs are interchangeable, providing ultimate flexibility.

A Base for Everything & Everything in its Place

The included multi-purpose Tono-Vera Vet Base is the perfect solution for docking and charging* your device while also conveniently storing and dispensing $Ocu-Dot^{\mathsf{T}}$ Tonometer Probes.

Ready for the Road

When your patients can't come to you, bring Tono-Vera Vet to them in the durable, lightweight, and protective Carry Case.

Electronic Health Record Ready

Tono-Vera Vet's unique Bluetooth® capability allows for software updates and seamlessly exports final IOP results to a computer.

*Rechargeable model only.



Tono-Vera Vet is designed, engineered, and assembled in the USA. We at Reichert® stand behind our products by backing up Tono-Vera™ Vet Tonometer with a two-year warranty.

Catalog Numbers

16309 Tono-Vera® Vet Tonometer Starter Kit - Rechargeable
16308 Tono-Vera® Vet Tonometer Starter Kit - AA Battery
16307 Ocu-Dot™ Tonometer Probes - 100 Count

Specifications

Dimensions: $19.8 \times 6.6 \times 3.0$ cm $(7.8 \times 2.6 \times 1.2$ in) Display Dimensions: 2.8×2.8 cm $(1.1 \times 1.1$ in) Weight (without batteries): 113.4g (4.0 oz)

Measurement Range: 1-99 mmHg

Measurement Distance Range: 6mm (±2mm)

Power Source: Lithium-ion or Four AA Batteries

Species: Dog, Cat, Horse, Rabbit **Languages:** English, German, French,

Spanish, Italian, Portuguese





Call Reichert toll-free at 1-888-849-8955 or +1 716-686-4500, or contact your Authorized Reichert Technologies® Distributor. 3362 Walden Avenue, Suite 100, Depew, NY 14043 USA ©2020 AMETEK, Inc. & Reichert Inc. 16308-110-

Reichert and Reichert Technologies are registered trademarks of Reichert, Inc. Tono-Vera, ActiVew, and Ocu-Dot are trademarks of Reichert, Inc. AMETEK is a registered trademark of AMETEK, Inc.

Bluetooth trademark is owned by the Bluetooth SIG.

passionate about eye care

